


**Configuración Cytoflex**

Configuración óptica  
**Cytoflex**  
 LNC-PNG03/R05-012


**1. LASER AZUL**



CANAL	FITC	PE	PC5.5	PC7
FLUOROCROMOS COMUNES	FITC	PE	PE Cy5.5	PE Cy7
	AF 488	PI	PE Cy5	PE AF750
	CFSE		PerCp	
	GFP		PerCp Cy5.5	


**Configuración Cytoflex**

**2. LASER ROJO**



CANAL	APC	APC A700	APC A750
FLUOROCROMOS COMUNES	APC	APC AF700	APC AF750
	AF 647	AF 700	APC Cy7
	EFluor 660		

**3. LASER VIOLETA**



CANAL	PB450	KO525	BV610
FLUOROCROMOS COMUNES	BV421	KO 525	BV 610
	PACIFIC BLUE	KO	BV 605
	PB450	BV510	Odot 605
		AmCyan	